

Preparation and Standardisation of Potassium Permanganate

Experiment 8 - Practical Notes

Name: _____ Roll no.: _____ Date: _____

Aim

To prepare a nominal titrant and determine its actual strength.

Principle

Potassium permanganate is an oxidant and a self-indicator. In acidic oxalate standardisation, two permanganate ions react with five oxalate ions; one mole permanganate accepts five electrons.

Requirements

- Assigned potassium permanganate solution and primary-standard sodium oxalate
- Prescribed sulphuric-acid medium, controlled water bath and burette

Procedure

1. Prepare

Prepare the nominal permanganate solution and complete the assigned conditioning and clarification before standardisation.

2. Prepare standard

Weigh primary-standard oxalate and dissolve in the prescribed sulphuric-acid medium.

3. Control conditions

Use the method-specified warming; the oxalate reaction is slow at the beginning.

4. Titrate

Add titrant at the prescribed rate and determine the persistent pale-pink endpoint.

5. Repeat

Obtain concordant titres, calculate strength and label the protected storage bottle.

Observations and Record

Reading	Trial 1	Trial 2	Trial 3
Standard mass, g	_____	_____	_____
Initial burette, mL	_____	_____	_____
Final burette, mL	_____	_____	_____
Corrected titre, mL	_____	_____	_____

Worked Example

The numbers below are illustrative; enter your own observations in the record.

- $N(KMnO_4) = \text{mass sodium oxalate (g)} / [67.00 \text{ g/equivalent} \times \text{titre (L)}]$.
- Example: 0.16750 g sodium oxalate and 25.00 mL gives 0.1000 N, equivalent to 0.0200 M permanganate in this acidic reaction.

Result

Record the observed outcome and calculated value, with units. State whether the stated laboratory specification was met; do not copy the worked example as your result.

Precautions

- Use the specified acid; hydrochloric acid can introduce side reactions.
- Do not use organic filter paper in an improvised clarification process.
- Match temperature and endpoint conditions across repeats.

Viva Questions and Answers

1. Why is $KMnO_4$ a self-indicator?

A slight excess imparts the endpoint colour.

2. Why is 0.1 N not 0.1 M here?

Each permanganate ion accepts five electrons in the acidic reaction.

3. Why use concordant readings?

Agreement reduces random endpoint and reading error, but does not remove systematic errors.

Continue Learning with FirstHope

Complete Notes on the FirstHope App

Complete subject and semesters are available on FirstHope app only.

- Get complete subject and semester notes in the FirstHope app.
- Available for B.Pharm, D.Pharm, M.Pharm and Pharm.D students.

[Click here to get the FirstHope app](#)

Join the FirstHope Telegram Channel

Join FirstHope Notes and keep up with channel updates.

[Click here to join Telegram](#)

Instagram

@firsthope.official

[Click here to follow Instagram](#)

YouTube Channel

FirstHope-Pharmacy / Live streams

[Click here to watch the YouTube channel](#)

WhatsApp Channel

Follow FirstHope on WhatsApp.

[Click here to join the WhatsApp channel](#)

FirstHope Website

firsthopenotes.com

[Click here to visit the website](#)