

Qualification of manufacturing equipment

Qualification of Manufacturing Equipment: Introduction

- **Qualification of manufacturing equipment** is a documented process that demonstrates equipment is properly designed, installed, operated, and capable of consistently performing according to predetermined specifications.
- It is a regulatory requirement under **Good Manufacturing Practices (GMP)** and forms the basis for process validation.

Objectives of Equipment Qualification

- Ensure equipment is suitable for its intended use.
- Confirm correct installation and operation.
- Verify consistent performance.
- Ensure compliance with GMP and regulatory guidelines.
- Produce safe and high-quality products.

Stages of Equipment Qualification

All manufacturing equipment is qualified in four stages:

1. Design Qualification (DQ)

- Verifies that the equipment design meets user requirements and GMP standards.
- Conducted before purchasing or installing equipment.
- **Activities:**
 - Review equipment specifications.
 - Check material of construction.
 - Verify production capacity.
 - Assess safety features.

2. Installation Qualification (IQ)

- Confirms that the equipment has been installed correctly according to manufacturer recommendations and approved drawings.
- **Activities:**
 - Verify equipment identification.
 - Check utilities (electricity, water, air, steam).
 - Confirm installation location.
 - Verify calibration status.
 - Review manuals and certificates.

3. Operational Qualification (OQ)

- Demonstrates that the equipment operates properly throughout its specified operating range.
- **Activities:**
 - Test alarms and controls.
 - Verify operating parameters.
 - Check safety devices.
 - Confirm operating limits.

4. Performance Qualification (PQ)

- Confirms that the equipment performs consistently under actual production conditions.
- **Activities:**
 - Manufacture trial batches.
 - Monitor product quality.
 - Verify reproducibility over multiple runs.

Qualification of Dry Powder Mixers

Purpose

- Used to mix powders uniformly before tablet or capsule manufacturing.

Critical Parameters

- Mixing time
- Mixing speed
- Blender load capacity
- Blend uniformity
- Mixing efficiency

Qualification

DQ (Design Qualification)

- Select appropriate mixer.
- Verify SS 316L construction.
- Confirm required batch size.

IQ (Installation Qualification)

- Verify installation and alignment.

- Check electrical supply.
- Confirm motor specifications.
- Verify timer and speed calibration.

OQ (Operational Qualification)

- Check rotation speed.
- Verify timer accuracy.
- Test emergency stop.
- Check operating controls.

PQ (Performance Qualification)

- Prepare three consecutive batches.
- Analyze samples from different mixer locations.
- Confirm uniform blend.

Acceptance Criteria

- Uniform blend achieved.
- No segregation.
- Equipment operates without malfunction.

Qualification of Fluid Bed Dryer (FBD)

Purpose

- Removes moisture from granules before tablet compression.

Critical Parameters

- Inlet air temperature
- Outlet air temperature
- Airflow rate
- Drying time
- Product moisture content

Qualification

DQ

- Verify drying capacity.
- Confirm heating system.

- Ensure GMP-compliant design.

IQ

- Check installation.
- Verify HEPA filters.
- Confirm utility connections.

OQ

- Test airflow.
- Verify temperature controls.
- Check alarms and safety interlocks.

PQ

- Dry three consecutive batches.
- Measure final moisture content.
- Ensure uniform drying.

Acceptance Criteria

- Uniform moisture content.
- Stable operating temperature.
- Consistent drying time.

Qualification of Tray Dryer

Purpose

- Removes moisture from granules by heated air circulation.

Critical Parameters

- Temperature distribution
- Air circulation
- Drying time
- Moisture removal

Qualification

- Verify temperature uniformity.
- Confirm airflow distribution.
- Monitor drying efficiency.
- Check thermostat accuracy.

Acceptance Criteria

- Uniform drying.
- Temperature within specified limits.

Qualification of Tablet Compression Machine

Purpose

- Compresses granules into tablets of required weight, hardness, thickness, and appearance.

Critical Parameters

- Compression force
- Turret speed
- Tablet weight
- Tablet thickness
- Tablet hardness

Qualification

DQ

- Verify production capacity.
- Confirm punches and dies.
- Review machine specifications.

IQ

- Check installation.
- Verify electrical supply.
- Confirm lubrication system.
- Review documentation.

OQ

- Test turret speed.
- Verify compression force.
- Check weight control system.
- Test emergency stop.

PQ

- Produce three validation batches and evaluate:



FREE PREVIEW

[Click here to visit our website](#)

PHARMACY PREMIUM NOTES

Complete subject
and semesters are
available on
FirstHope app only.



[Scan or click here](#)

Click here to get the app



Available on Google Play

JOIN THE FIRSTHOPE COMMUNITY

Join our Telegram channel

Click here to join Telegram



[Scan or click here](#)

@FirstHopeNotes

Instagram

@firsthope.official

[Click here to follow >](#)

YouTube channel

FirstHope-Pharmacy / Live streams

[Click here to watch >](#)

WhatsApp channel

Follow FirstHope on WhatsApp

[Click here to join >](#)

FirstHope website

firsthopenotes.com

[Click here to visit >](#)

B.Pharm / D.Pharm / M.Pharm / Pharm.D

Preview notes. Complete subjects and semesters in the FirstHope app.