

High-Performance Thin Layer Chromatography (HPTLC)

High-Performance Thin Layer Chromatography (HPTLC)

Introduction

- High-Performance Thin Layer Chromatography (HPTLC) is an advanced and automated form of traditional Thin Layer Chromatography (TLC).
- It offers greater resolution, sensitivity, reproducibility, and speed, making it suitable for qualitative and quantitative analysis in pharmaceutical and biomedical research.
- HPTLC is widely accepted by pharmacopeias (USP, IP, BP) and regulatory authorities for:
 - Identification of active pharmaceutical ingredients (APIs)
 - Estimation of impurities and degradation products
 - Quality control of herbal formulations
- It is especially useful in multi-component analysis and routine quality control where high sample throughput is required.

Theory of HPTLC

- HPTLC is based on the same fundamental principle as TLC:
 - Partition and adsorption chromatography, where separation occurs due to differences in the affinities of components toward the stationary and mobile phases.
- **In HPTLC:**
 - The stationary phase is a high-resolution, precoated layer of fine particle silica or alumina (5–10 μm) with a uniform and dense coating.
 - The mobile phase is a solvent or solvent mixture that moves via capillary action.
 - The components of a sample partition between the two phases and separate based on their polarity, solubility, and molecular interaction.
- Because of its controlled environment and automation, HPTLC allows precise migration, well-defined R_f values, and quantitative data output.

Instrumentation of HPTLC

The HPTLC system consists of modular, computer-controlled units, which offer high precision and reproducibility.

1. Sample Applicator (Automatic Spotter)

- Precisely applies the sample as narrow bands or spots using a micro-syringe or capillary system.
- Ensures uniformity, which is critical for quantitative analysis.
- Multiple samples can be spotted in parallel on the same plate.

2. Precoated HPTLC Plates

- Glass, aluminum, or plastic-backed plates with a uniform layer of adsorbent (e.g., silica gel 60 F254).
- Particle size: 5–10 μm → ensures high resolution.
- **Plate size:** typically 10 × 10 cm or 20 × 10 cm.

3. Development Chamber (Automated or Manual)

- Saturated with mobile phase vapor before development.
- Ensures uniform solvent migration.
- Can use isocratic or gradient solvent systems.

4. Densitometer/Scanner

- Scans developed chromatograms at specific wavelengths (UV/Visible).
- Measures peak area or height → directly proportional to concentration.
- Enables quantification of separated drugs.

5. Data Station / Software

- Controls all modules: sample application, plate development, scanning.
- Performs integration, calibration, and quantitative analysis.
- Stores data for regulatory compliance and audit trails.

6. Derivatization Apparatus (Optional)

- Used to visualize colorless compounds via post-development chemical reactions.
- Example: Ninhydrin for amino acids, Dragendorff's reagent for alkaloids.

Applications of HPTLC in Pharmaceutical Analysis

A. Qualitative Analysis

- Identification of drugs by comparing R_f and spectra with standards.
- **Example:** Detecting aspirin in combined pain-relief tablets.

B. Quantitative Analysis

- Measures concentration of drug substances or impurities using peak area/height.
- Suitable for bulk drug analysis, formulations, and dissolution testing.

C. Impurity Profiling

- Detects related substances, residual solvents, and degradation products.
- Stability-indicating method for routine testing.
- **Example:** Monitoring photodegradation of nifedipine under UV exposure.

D. Drug-Excipient Compatibility

- Can determine if excipients affect the stability or recovery of drugs.
- **Example:** Separation of ibuprofen from lactose, magnesium stearate, and talc.

E. Herbal Drug Standardization

- Simultaneous identification and quantification of multiple phytoconstituents.
- Fingerprinting helps in authenticity, purity, and quality control.
- **Example:** Quantification of curcumin in turmeric extract.

F. Dissolution and Content Uniformity Testing

- Used in tablets and capsules to confirm uniform distribution of the drug.
- Validates consistency in manufacturing.

G. Bioanalytical Applications

- Though not as sensitive as LC-MS or HPLC, HPTLC can analyze biological fluids for drugs and metabolites.
- **Example:** Screening urine samples for paracetamol or caffeine.

H. Cost-effective Screening Tool

- Ideal for screening counterfeit or substandard drugs, especially in low-resource settings.
- Can test multiple samples simultaneously, making it efficient for QC labs.



FREE PREVIEW

[Click here to visit our website](#)

PHARMACY PREMIUM NOTES

Complete subject
and semesters are
available on
FirstHope app only.



[Scan or click here](#)

Click here to get the app



Available on Google Play

JOIN THE FIRSTHOPE COMMUNITY

Join our Telegram channel

Click here to join Telegram



[Scan or click here](#)

@FirstHopeNotes

Instagram

@firsthope.official

[Click here to follow >](#)

YouTube channel

FirstHope-Pharmacy / Live streams

[Click here to watch >](#)

WhatsApp channel

Follow FirstHope on WhatsApp

[Click here to join >](#)

FirstHope website

firsthopenotes.com

[Click here to visit >](#)

B.Pharm / D.Pharm / M.Pharm / Pharm.D

Preview notes. Complete subjects and semesters in the FirstHope app.