

Separation of Sugars by Thin-Layer Chromatography

Practical 09 | Instrumental Methods of Analysis and Novel Drug Delivery Systems

Aim

To study Separation of Sugars by Thin-Layer Chromatography, and interpret the recorded observations.

Principle

- The assigned TLC system separates sugar standards and a mixture for migration-pattern comparison. Detection and R_f depend on the stationary/mobile-phase and visualization system.
- Compounds separate because they interact differently with the stationary and mobile phases.
- Identity assessment uses an appropriate comparison standard and method conditions; one migration or retention value alone is not conclusive.

Materials and Apparatus

- Prepared sugar standards/mixture, assigned TLC plate and development/detection system.
- Assigned test extract or mixture, comparison standards and the specified stationary/mobile phase system.
- Appropriate chromatographic apparatus, application equipment, developing chamber or instrument and the assigned detector.

Method

1. Record the stationary phase, mobile phase and the purpose of the separation.
2. Prepare test and comparison solutions at suitable concentrations for the assigned method.
3. Introduce measured test and comparison portions without overloading the separation system.
4. Develop or run the separation under the assigned conditions, recording the relevant operating details.
5. Locate separated bands, spots or peaks using the specified observation method.
6. Compare the test pattern with the standard and calculate the appropriate separation metric.

Observation Record

Complete this record using actual laboratory findings or the provided teaching data.

Observation / parameter	Record / finding
Sugar/mixture	_____
Spot appearance	_____
Spot distance	_____
Front distance	_____
Rf and match	_____

Calculations and Interpretation

- For planar chromatography: $R_f = \text{distance travelled by spot centre} / \text{distance travelled by solvent front}$.
- R_f is dimensionless and normally lies between zero and one under a valid development.

Worked Illustration

Illustration: a spot travelling 3.0 cm with a solvent front travelling 6.0 cm has $R_f = 0.50$.

The values in this illustration are for calculation practice.

Result

Record the actual finding, calculated value with units where relevant, and a conclusion supported by the observations.

Finding: _____

Conclusion: _____

Precautions

- Keep application concentrated and avoid overloading.
- Record the solvent front immediately for planar chromatography.
- Use retention or migration data only with the stated system and relevant identity evidence.

Viva Questions

1. Why does overloading impair separation?

It broadens bands or spots and can change apparent migration or retention.

2. Why run a comparison standard?

It supplies a relevant reference pattern under the same analytical conditions.

Continue Learning with FirstHope

Complete Notes on the FirstHope App

Complete subject and semesters are available on FirstHope app only.

- Get complete subject and semester notes in the FirstHope app.
- Available for B.Pharm, D.Pharm, M.Pharm and Pharm.D students.

[Click here to get the FirstHope app](#)

Join the FirstHope Telegram Channel

Join FirstHope Notes and keep up with channel updates.

[Click here to join Telegram](#)

Instagram

@firsthope.official

[Click here to follow Instagram](#)

YouTube Channel

FirstHope-Pharmacy / Live streams

[Click here to watch the YouTube channel](#)

WhatsApp Channel

Follow FirstHope on WhatsApp.

[Click here to join the WhatsApp channel](#)

FirstHope Website

firsthopenotes.com

[Click here to visit the website](#)